

# BTI EXTRACTORS COMPATIBILITY CHART

IMPLANT BRAND	MODEL	✓	Nº
<b>BTI</b>	INTERNA	✓	<b>3</b>
	EXTERNA	✓	<b>1</b>
	TINY	✓	<b>1</b>
<b>BIOMET 3i</b>	NT Certain <sup>®</sup>	✓	<b>1</b>
	Certain <sup>®</sup> /Prevail <sup>®</sup>	✓	<b>3</b>
	Hexágono Externo	✓	<b>1</b>
<b>ASTRA TECH</b>	SMALL (Ø3,5/4)	✓	<b>0</b>
	LARGE (Ø4,5/5)	✓	<b>3</b>
<b>STRAUMANN</b>	Regular Neck (RN)	✓	<b>3</b>
<b>NOBEL BIOCARE</b>	Replace Ø3,5	✓	<b>2</b>
	Replace Ø4,3	✓	<b>2 · 3</b>
	Replace Ø5	✓	<b>4</b>
	Branemark System NP	✓	<b>1</b>
	Branemark System RP	✓	<b>1</b>
<b>BICON</b>	Narrow Implants Ø 3,5	✓	<b>1</b>
	Narrow Implants Ø 4	✓	<b>1</b>
	Short implants Ø 4,5	✓	<b>3</b>
	Short implants Ø 5	✓	<b>3</b>
	Short implants Ø 6	✓	<b>3</b>
<b>BEGO SEMADOS</b>	S Line Ø 3,25	✓	<b>1 · 2</b>
<b>DENTSPLY FRIADENT</b>	XIVE <sup>®</sup> S D3	✓	<b>1 · 2</b>
	XIVE <sup>®</sup> S D3,8	✓	<b>2</b>
	XIVE <sup>®</sup> S D5,5	✓	<b>4</b>
	ANKYLOS <sup>®</sup> C/X A Ø 3,5	✓	<b>1</b>
	ANKYLOS <sup>®</sup> C/X B Ø 4,5	✓	<b>1</b>
	ANKYLOS <sup>®</sup> C/X C Ø 5,5	✓	<b>1</b>
	ANKYLOS <sup>®</sup> C/X D Ø 7	✓	<b>1</b>
	FRIALIT <sup>®</sup> D3,8	✓	<b>2</b>
	FRIALIT <sup>®</sup> D4,5	✓	<b>2</b>
	FRIALIT <sup>®</sup> D5,5	✓	<b>2</b>
	FRIALIT <sup>®</sup> D6,5	✓	<b>2</b>
<b>LIFECORE</b>	RENOVA <sup>™</sup> SDI	✓	<b>1</b>
	RENOVA <sup>™</sup> RDI	✓	<b>1</b>
<b>NEODENT</b>	CM	✓	<b>1</b>
	GT	✓	<b>1</b>
	Titamax Poros Ø 3,3	✓	<b>1</b>
	Titamax Alvim Ø 3,5	✓	<b>1</b>
	Titamax Ti Ø 3,75 – 4,0	✓	<b>1</b>
Alvim Ti Ø 4,3	✓	<b>1</b>	
<b>MIS</b>	SEVEN 4,2	✓	<b>2 · 3</b>
<b>KLOCKNER</b>	S3M 3112	✓	<b>0</b>
	Cone Essential	✓	<b>3</b>
	Cone 1.5 Essential	✓	<b>3</b>

Nº 0 · EXIME

Max. 70 Ncm of torque



Nº 1 · EXIM



Nº 2 · INEXIML



Nº 3 · INEXIM



Nº 4 · EXIMA



✓ OPTIMUM COMPATIBILITY TESTED BY BTI

# BTI EXTRACTORS COMPATIBILITY CHART

IMPLANT BRAND	MODEL	✓	Nº
<b>MOZOGRAU</b>	MG-OSSEOUS STD Ø 3,4	✓	<b>1</b>
	MG-OSSEOUS STD Ø 3,75	✓	<b>1</b>
	MG-OSSEOUS STD Ø 4,25	✓	<b>1</b>
<b>IMPLANT DIRECT LLC</b>	Screw Direct Ø 3,7	✓	<b>1</b>
	Screw Direct Ø 4,7	✓	<b>1</b>
<b>CAMLOG</b>	CAMLOG ROOT LINE Ø 5	✓	<b>4</b>
	CAMLOG SCREW LINE Ø 5	✓	<b>4</b>
	CONELOG SCREW LINE Ø 5	✓	<b>4</b>
<b>SOUTHERN IMPLANTS</b>	BAT	✓	<b>1</b>
<b>PHIBO</b>	TSA ADVANCE S3	✓	<b>0</b>
	TSA ADVANCE S5	✓	<b>1</b>
	TSA S3	✓	<b>1</b>
	TSA S4	✓	<b>1</b>
	TSH S3	✓	<b>1</b>
	TSH S4	✓	<b>1</b>
	BNT S2	✓	<b>1*</b>
<b>DENTAURUM IMPLANTS</b>	TioLogic	✓	<b>0 - 1</b>

(\*) This 3.3 diameter implant can not resist the torque applied by the LLT200 wrench so the LLMQ surgical wrench should be used in this case.

This document shows confidential information and it can only be in BTI employees hands. It is absolutely forbidden to deliver, distribute, publish or any other misuse of it. If so it will be considered a serious breach of the non disclosure agreement and BTI will be entitled to take any legal actions in its hands to defend its rights.

## Nº 0 · EXIME

Max. 70 Ncm of torque



## Nº 1 · EXIM



## Nº 2 · INEXIML



## Nº 3 · INEXIM



## Nº 4 · EXIMA

